




Silicone	Coiled	15F
Material	Tip Design	French Size



COILED TIP

Provides a gentler, less traumatic, showering effect of dialysate on the peritoneum and prevents catheter migration.



LARGE SIDE HOLES

Increase flow and reduce incidence of omental occlusion.

TITANIUM BARBED CONNECTOR

Securely joins the two catheter end segments for upper abdominal or Presternal placement.



CATHETER ONLY		
MC20XC	15F X 64CM X-SERIES, COILED	5/BOX
BASIC SET		
MC20XS	15F X 64CM X-SERIES, COILED	5/BOX
ACCESSORIES		
MP217	17F VASCU-SHEATH INTRODUCER SET	5/BOX
TRX	X-SERIES TUNNELING ROD	5/BOX
MTC	TITANIUM CONNECTOR	5/BOX
FSX	X-SERIES FALLER STYLET	5/BOX
40501	PREOPERATIVE MARKING STENCIL PACK	1/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter
- (1) Titanium Barbed Connector
- (1) Female Luer
- (1) Clamp
- (1) Injection Cap

BASIC SET CONTAINS:

- (1) Catheter
- (1) Titanium Barbed Connector
- (1) Female Luer
- (1) Clamp
- (1) Injection Cap
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Peelable Sheath Introducer
- (1) Tunneler
- (1) 10cc Syringe
- (6) 4" x 4" Gauze
- (1) Injection Cap

